

FIG. 24

← Illuminating and Imaging Section

Excitation beam path

- L1: 50 mm f.l.
- L2: 150 mm f.l.
- L3: cylindrical lens 75 mm f.l.
- S1: silt aperture, 25 μ m wide x 3 mm long
- L4: 75 mm f.l.
- BS: beam splitter (dichroic if used in fluorescence)
- M1*: front side of bilateral scanning mirror
- L5: 150 mm f.l.
- W1 & W2: glass windows index matched to fiber bundle (not needed for fluorescence)
- OL: objective lens

Return beam path

- L6: 75 mm f.l.
- S2: silt aperture, 25 μ m wide x 3 mm long
- L7: 75 mm f.l.
- M1**: back side of bilateral scanning mirror
- L8: 150 mm f.l.

Approved
2/20/88
LD

Illuminating and Imaging Section ↙





FIG. 28

- Excitation beam path**
- L1: 50 mm f.l.
 - L2: 150 mm f.l.
 - L3: cylindrical lens 75 mm f.l.
 - S1: slit aperture, 25 μ m wide x 3 mm long
 - L4: 75 mm f.l.
 - BS: beam splitter (dichroic if used in fluorescence)
 - M1*: front side of bilateral scanning mirror
 - L5: 150 mm f.l.
 - W1 & W2: glass windows index matched to fiber bundle (not needed for fluorescence)
 - OL: objective lens
- Return beam path**
- L6: 75 mm f.l.
 - S2: slit aperture, 25 μ m wide x 3 mm long
 - DET: linear detector array or line scan camera



Probe Section

← Illuminating and Imaging Section

Excitation beam path

- L1: 50 mm f.l.
- L2: 150 mm f.l.
- L3: cylindrical lens 75 mm f.l.
- S1: slit aperture, 25 μ m wide x 3 mm long
- L4: 75 mm f.l.
- BS: beam splitter (dichroic II used in fluorescence)
- M1*: front side of bilateral scanning mirror
- L5: 150 mm f.l.
- W1 & W2: glass windows index matched to fiber bundle (not needed for fluorescence)
- FB1: fiber bundle
- OL: objective lens

Return beam path

- L6: 75 mm f.l.
- S2: slit aperture, 25 μ m wide x 3 mm long
- L7: 75 mm f.l.
- M1**: back side of bilateral scanning mirror
- L8: 150 mm f.l.
- FB2: identically mapped fiber bundle as FB1

FIG. 2C